


Clearline tips 1250 µL Hybridrack ClearLine® sterile low retention

- In clear polypropylene
- Store in a clean and dry place, at room temperature
- Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	258365
Packaging	HybridRack: cardboard rack with recyclable tip tray. (8 racks of 96 tips)
Volume	125 µl
Graduation	10, 50, 100, 200
Length	102 mm
Color	Natural
Filter	No
Low Retention	Yes
Sterility	Yes
Autoclavable	Yes



Units :
768 Tips