


**200 µl filter tip, LTS compatible, graduated, sterile, MultiRack, yellow tray ClearLine**

- Volume : 200 µl
- Longueur : 50,2 mm
- Packaging : In rack Clearline Plus - 10 x 96
- Compatibility : Compatible with Rainin LTS™ micropipettes
- Tray color : yellow
- Low retention : no
- Autoclavable : no
- Sterile : yes
- Filter : yes
- Color : natural
- Graduation : 10, 50, 100
- Material : polypropylene



**Units : 960**