


**1000 µl filter tip, LTS compatible, graduated, sterile, MultiRack, blue tray ClearLine**

- Volume : 1000 µl
- Longueur : 98,6 mm
- Packaging : In rack Clearline Plus - 8 x 96
- Compatibility : Compatible with Rainin LTS™ micropipettes
- Tray color : blue
- Low retention : no
- Autoclavable : no
- Sterile : yes
- Filter : yes
- Color : natural
- Graduation : 200, 500, 1000
- Material : polypropylene



**Units : 768**