


# Terumo® 3-part Syringe

## Product Information Card



### Description:

Sterile 3-part syringe for single use of *Terumo (Philippines) Corporation origin*.

The syringe is intended for the injection of fluids - immediately after filling - into the human body, for aspiration of fluids or for the sampling of blood or other body fluids.

### Components:

Description	Raw Materials
Barrel	Polypropylene
Plunger	Polystyrene (for 1 ml syringe only) Polypropylene
Gasket	Thermoplastic elastomer
Lubricant	Silicone oil



### Product Range:

Product Code	Size	Tip	Replacement for
SS+01T1	1 ml	Concentric, Luer Slip	BS-01T
SS+02S1	2.5 ml	Concentric, Luer Slip	SS*02SE1
SS+03L1	3 ml	Concentric, Luer Lock	
SS+05S1	5 ml	Concentric, Luer Slip	
SS+10ES1	10 ml	Eccentric, Luer Slip	SS-10ES
SS+20ES1	20 ml	Eccentric, Luer Slip	
SS+20L1	20 ml	Concentric, Luer Lock	
SS+50ES1	50 ml	Eccentric, Luer Slip	BS-50ES
SS+50L1	50 ml	Concentric, Luer Lock	BS-50LG
SS+50C1	50 ml	Catheter Tip	BS-50C

### Technical Data:

Product Code	Nominal Capacity	Graduated Capacity	Increments
SS+01T1	1 ml	1 ml	0.01
SS+02S1	2.5 ml	2.5 ml	0.1
SS+03L1	3 ml	3 ml	0.1
SS+05S1	5 ml	6 ml	0.2
SS+10ES1	10 ml	10 ml	0.2
SS+20ES1	20 ml	20 ml	1.0
SS+20L1	20 ml	20 ml	1.0
SS+50ES1	50 ml	60 ml	1.0
SS+50L1	50 ml	60 ml	1.0
SS+50C1	50 ml	60 ml	1.0

# Terumo® 3-Part Syringe

## Packaging:

*Single unit:* blister pack (film/film)

Product Code	Unit Box	Shipping Carton		
		Quantity	Gross Weight	Dimensions (LxWxH)
SS+01T1	100	1800	7.50 kg	430 x 418 x 306 mm
SS+02S1	100	2400	7.00 kg	496 x 427 x 427 mm
SS+03L1	100	1800	7.50 kg	640 x 385 x 262 mm
SS+05S1	100	1200	9.40 kg	480 x 379 x 362 mm
SS+10ES1	100	1200	11.00 kg	546 x 454 x 428 mm
SS+20ES1	50	600	10.15 kg	570 x 400 x 366 mm
SS+20L1	50	600	10.15 kg	570 x 400 x 366 mm
SS+50ES1	25	100	5.10 kg	466 x 312 x 268 mm
SS+50L1	25	100	5.10 kg	466 x 312 x 268 mm
SS+50C1	25	100	5.40 kg	466 x 312 x 268 mm

## Sterilization:

With e-beam radiation to a Sterility Assurance Level of at least 10<sup>-6</sup>

## Storage:

Shelf life of 5 years

\* The batch number and expiry date are embossed in the film of the blister pack

\* Do not store at extreme temperature and humidity

## Precautions for Use:

Check integrity of the unit package before use

For single use only

Dispose of safely after single use to avoid risk of infection

## Status of Regulation:

Class IIa, sterile product following MDD 93/42/EEC

Notified body:

TÜV Rheinland 0197

Manufacturer:

TERUMO (PHILIPPINES) CORPORATION

EC representative:

TERUMO EUROPE NV

Interleuvenlaan 40

B - 3001 Leuven

Belgium



**Terumo Corporation**  
44-1, 2-chome, Hatagaya  
Shibuya-ku, Tokyo 151-0072, Japan  
Tel.: +81 3 3374 8111  
Fax: +81 3 3374 8196

**Terumo Europe NV**  
Interleuvenlaan 40  
3001 Leuven, Belgium  
Tel.: +32 16 38 12 11  
Fax: +32 16 40 02 49

**Terumo Europe NV  
Benelux Sales Division**  
Belgium:  
Tel.: +32 16 39 25 80  
The Netherlands:  
Tel.: 0800 0220396

**Terumo Europe NV  
Emerging Markets Division**  
Tel.: +32 16 38 13 08

**Terumo Deutschland GmbH**  
Tel.: +49 6196 80 230

**Terumo Deutschland GmbH  
Zweigniederlassung  
Spreitenbach**  
Tel.: +41 56 419 10 10

**Laboratoires  
Terumo France S.A.S.**  
Tel.: +33 130 96 13 00

**Terumo Italia S.r.l.**  
Tel.: +39 0651 96 14 20

**Terumo Europe España SL**  
Tel.: +34 902 10 12 98

**Terumo UK Ltd**  
Tel.: +44 1276 480 440

**Terumo Russia LLC**  
Tel.: + 7 495 988 4740

**Terumo Sweden AB**  
Tel.: +46 3174 85 880

**Terumo Denmark**  
Tel.: +45 7020 93 80

**Terumo Corporation Dubai Branch  
(Middle East and Africa)**  
Tel.: +971 4 221 2220

**Terumo BCT Tıbbi Cihazlar Dağıtım  
ve Hizmetleri A.Ş.**  
Tel.: +90 216 590 0074



HO 20GB-0307TE-VIII(01.15)E

© Registered Trademark  
Published by Terumo Europe NV