



TECHNICAL DATASHEETPART NUMBER13245Image: Data of the second s

ClearLine Tips

- In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	713245
Packaging	ClearLine plus rack with hinged lid
Volume	10 µl
Graduation	1, 2, 5, 10
Length	44,8 mm
Color	Natural
Filter	No
Low Retention	No
Sterility	No
Autoclavable	Yes

Sales units :

960 tips

Dominique Dutscher SAS 30 rue de l'Industrie ZAC BRUMATH NORD BP62 67172 BRUMATH CEDEX Tél : +33 (0)3 88 59 33 90 Fax : +33 (0)3 88 59 33 99 Mail : info@dutscher.com Web : www.dutscher.com