



TECHNICAL DATASHEET

PART NUMBER

713250





The ISO8655 norm guarantees that the tips used do not alter the technical characterics of the pipettes, as given by the manufacturer.










ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified: RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability: 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	713250
Packaging	ClearLine plus rack with hinged lid
Volume	200 μΙ
Graduation	10, 50, 100
Length	51,44 mm
Color	Natural
Filter	No
Low Retention	Yes
Sterility	Yes
Autoclavable	Yes
Filter Low Retention Sterility	No Yes Yes



Sales units:

960 tips