


## TECHNICAL DATASHEET

## PART NUMBER

**713265**




The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.



### ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	713265
Packaging	ClearLine plus rack with hinged lid
Volume	300 µl
Graduation	10, 50, 100, 200
Length	77,5 mm
Color	Natural
Filter	Yes
Low Retention	Yes
Sterility	Yes
Autoclavable	No



### Sales units :

960 tips