






Designed for reliable labelling in laboratories

edding 8014 / 8015 F laboratory marker

- Suitable for fine, precise labelling of smooth and transparent surfaces, such as slides, reaction vessels, petri and laboratory dishes
- Recommended by Clariant International AG
- Extremely quick-drying, wear-resistant and waterproof, lightfast and low-odour ink
- Permanent on plastics, removable from glass with acetone, ethanol, isopropanol and methyl ethyl ketone
- Excellent frost (until $-183\text{ }^{\circ}\text{C}$), heat (on aluminium until $+500\text{ }^{\circ}\text{C}$) and abrasion resistant (see TÜV certificate)



Order Info p. 14



Order Info p. 14

* Certification No. 2716 for edding 8014 laboratory marker / Certification No. 2736 for edding 8015 F laboratory marker