


CAT. NO.....:	162-0020/-0045		161-0020/-0045	
PORE SIZE (µm).....:	0.2/0.45		0.2/0.45	
UPPER CAP. (mL).....:	500		1000	
RECEIVER CAP. (mL)::	500		1000	
MEMBRANE COLOR.....:	WHITE		WHITE	
MEMBRANE DIA (mm)::	90		90	
	INCHES	MM	INCHES	MM
A (UPPER I.D.).....:	3.84	97.5	4.48	113.8
B HEIGHT.....:	10.7	272	12.8	325
C DIA (BASE).....:	3.86	98.0	4.48	113.8
D OVERALL.....:	4.0	102	4.2	106
NOM. WT. (gms).....:	170		250	
MATERIALS				
UPPER CONTAINER	PS			
COVER	PS			
COLLAR	YELLOW PS			
RECEIVER	PS			
STORAGE CLOSURE	WHITE HDPE			
ADAPTOR	WHITE LDPE			
GASKET	FOAMED LDPE			
MEMBRANE	SFCA			
VENT PLUG	CA			



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Ponorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MG/MC	APP'D	PVC	PSS 00-016X-XXXX-2	REV.	0699
CK'D	PVC		6/25/99			

DESCRIPTION: SERUM AND MEDIA-PLUS  
SFCA FILTER UNITS