


CAT. NO.....:	165-0045	168-0045	166-0045			
PORE SIZE (μm).....:	0.45	0.45	0.45			
UPPER CAP. (mL).....:	150	250	500			
RECEIVER CAP. (mL).....:	150	250	500			
MEMBRANE COLOR.....:	WHITE	WHITE	WHITE			
MEMBRANE DIA (mm).....:	50	50	75			
	INCHES	MM	INCHES	MM	INCHES	MM
A (UPPER I.D.).....:	3.15	80.0	3.15	80.0	3.80	96.5
B HEIGHT.....:	6.15	156.2	7.9	200	9.8	249
C DIA (BASE).....:	3.50	88.9	3.50	88.9	3.86	98.0
D OVERALL.....:	3.7	94	3.7	94	4.0	102
NOM. WT. (gms).....:	78	95	155			

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	BLUE PS
RECEIVER	PS
STORAGE CLOSURE	WHITE HDPE
ADAPTOR	WHITE LDPE
GASKET	FOAMED LDPE
MEMBRANE	PES
VENT PLUG	CA

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MG/MC	APP'D	PVC	PSS 00-016X-0045-1	REV.	1205
CK'D	SH	12/22/05				

DESCRIPTION:

PES FILTER UNITS