

Chest freezer XLT C75



- Low cost of ownership
- Limited maintenance and on-site service
- Easily installable & movable

- Efficient cyclopentane insulation
- Eco-friendly, HCFC-free Nature R refrigerant
- ISO 9001, 14001 & 45001 Certified

Specifications

Technical

| | |
|--------------------------|------------------|
| Gross volume | 74 litres |
| Net Volume | 71 litres |
| Temp. range | -40°C to -60°C |
| Temp. probes | 1/PT 1000 |
| Compressor | Embraco |
| Baskets (incl./optional) | 1/2 |
| Refrigerant | Nature R |
| Noise level | 55 dB |
| Power consumption 24/h | 4,0 kWh |
| Power consumption | 1.460 kWh / year |
| Voltage | 230V |
| Frequency | 50 Hz |

Features

| | |
|--|----------|
| Audible and visual alarm | YES |
| Castors | NO |
| Integrated lock | YES |
| Inner lids | NO |
| Defrost | Manual |
| Battery Backup | NO |
| USB | NO |
| Data connection | NO |
| Digital controller with display | XR30CX |
| Data logging | NO |
| High/Low temperature alarm | YES |
| Porthole for external probe & CO2 backup | Ø12 mm |
| Connection for remote alarm | Optional |

Material

| | |
|-------|-----------------------------|
| Inner | Precoated white steel |
| Outer | Galvanized pre-coated steel |

Outer dimensions

| | |
|----------------------------------|--------|
| Height | 870 mm |
| Width | 555 mm |
| Depth (incl. condenser) | 580 mm |
| Depth (incl. condenser & handle) | 630 mm |
| Net weight | 48 kg |

Inner dimensions

| | |
|-------------------|--------|
| Height | 466 mm |
| Width | 387 mm |
| Depth | 387 mm |
| Insulation | 80 mm |
| Insulation type | PUR |
| Internal humidity | <5 % |

XLT C75



Product no.

N112001

Storage capacity

| Box height | No. of boxes | No. of racks | Total no. of boxes | Product no. |
|------------|--------------|--------------|--------------------|-------------|
| 50 mm | 8 | 5 | 40 | N311210 |
| 75 mm | 5 | 5 | 25 | N311310 |
| 100 mm | 4 | 5 | 20 | N311410 |