

COOLSAFE TOUCH 55-4

SPECIFICATIONS



For more information please visit:
www.labogene.com

CoolSafe Touch 55-4	Unit	Part no: 700047
Digital read-out of temperature		Yes
Digital read-out of pressure		Yes
Ultimate temperature at 20°C room temperature	°C	-55
Power requirements	V/Hz	230/50-60
Power consumption	W	800
Insulation	mm	90
Condenser capacity	kg	3
Condenser capacity per 24 hours	kg	2.5
Total volume of Condenser	L	4
Weight net	kg	40
Cabinet dimensions (DxWxH)	mm	500x400x520
Cabinet material		Polyester coated steel
Condenser material		Stainless steel AISI 316
Cooling media		R 507
Vacuum Valve		Electric
Drain tap		Yes
Touch screen for control options		Yes
USB port for LaboGene software/SD card for data		Yes

LaboGene ApS
 Bjarkesvej 5
 DK-3450 Allerød

Tel (+45) 3940 2566
 Fax (+45) 3995 2566
 Mail info@labogene.com
 Web www.labogene.com

