

COOLSAFE TOUCH 110-4

SPECIFICATIONS



For more information please visit:
www.labogene.com

| CoolSafe Touch 110-4 | Unit | Part no: 700049 |
|---|------|--------------------------|
| Digital read-out of temperature | | Yes |
| Digital read-out of pressure | | Yes |
| Ultimate temperature at 20°C room temperature | °C | -110 |
| Power requirements | V/Hz | 230/50-60 |
| Power consumption | W | 1300 |
| Insulation | mm | 90 |
| Condenser capacity | kg | 3 |
| Condenser capacity per 24 hours | kg | 2.5 |
| Total volume of Condenser | L | 4 |
| Weight net | kg | 55 |
| Cabinet dimensions (DxWxH) | mm | 500x400x520 |
| Cabinet material | | Polyester coated steel |
| Condenser material | | Stainless steel AISI 316 |
| Cooling media | | R 507/R1150 |
| Vacuum Valve | | Electric |
| Drain tap | | Yes |
| Touch screen for control options | | Yes |
| USB port for LaboGene software/SD card for data | | Yes |

LaboGene ApS
 Bjarkesvej 5
 DK-3450 Allerød

Tel (+45) 3940 2566
 Fax (+45) 3995 2566
 Mail info@labogene.com
 Web www.labogene.com

