

# COOLSAFE SUPERIOR TOUCH XS SPECIFICATIONS



For more information please visit:  
[www.labogene.com](http://www.labogene.com)

CoolSafe Superior Touch XS	Unit	Part no: 700059
Digital read-out of temperature		Yes
Digital read-out of pressure		Yes
Ultimate temperature at 20°C room temperature	°C	-100
Power requirements	V/Hz	230/50-60
Power consumption	W	1300
Insulation	mm	90
Condenser capacity when equipped with the rack	kg	1
Condenser capacity per 24 hours when equipped with the rack	kg	1
Total volume of Condenser	L	9
Weight net	kg	75
Cabinet dimensions (DxWxH)	mm	620x495x745
Cabinet material		Polyester coated steel
Condenser material		Stainless steel AISI 316
Cooling media		R 507/R170
Vacuum Valve		Electric
Drain tap		Yes
Touch screen for control options		Yes
USB port for LaboGene software/SD card for data		Yes
3 electric heated shelves ø180, for pre freeze and freeze drying		Yes, 0,25L/shelf
Distance between shelves	mm	34

LaboGene ApS  
Bjarkesvej 5  
DK-3450 Allerød

Tel (+45) 3940 2566  
Fax (+45) 3995 2566  
Mail [info@labogene.com](mailto:info@labogene.com)  
Web [www.labogene.com](http://www.labogene.com)

