

COOLSAFE SUPERIOR TOUCH XL 3 SHELVES SPECIFICATIONS



For more information please visit:
www.labogene.com

CoolSafe Superior Touch XL 3 Shelves	Unit	Part no: 700061
Digital read-out of temperature		Yes
Digital read-out of pressure		Yes
Ultimate temperature at 20°C room temperature	°C	-95
Power requirements	V/Hz	230/50-60
Power consumption	W	1400
Insulation	mm	90
Condenser capacity when equipped with the rack	kg	2
Condenser capacity per 24 hours when equipped with the rack	kg	2
Total volume of Condenser	L	15
Weight net	kg	85
Cabinet dimensions (DxWxH)	mm	620x495x745
Cabinet material		Polyester coated steel
Condenser material		Stainless steel AISI 316
Cooling media		R 507/R170
Vacuum Valve		Electric
Drain tap		Yes
Touch screen for control options		Yes
USB port for LaboGene software/SD card for data		Yes
3 electric heated shelves ø180, for pre freeze and freeze drying		Yes, 0,25L/shelf
Distance between shelves	mm	79

LaboGene ApS
 Bjarkesvej 5
 DK-3450 Allerød

Tel (+45) 3940 2566
 Fax (+45) 3995 2566
 Mail info@labogene.com
 Web www.labogene.com

