

COOLSAFE PRO 110-4

SPECIFICATIONS



For more information please visit:
www.labogene.com

CoolSafe Pro 110-4	Unit	Part no: 7001000516
Digital read-out of temperature		Yes
Digital read-out of pressure		Yes
Ultimate temperature at 20°C room temperature	°C	-110
Power requirements	V/Hz	230/50-60
Power consumption	W	600
Insulation	mm	90
Condenser capacity	kg	3
Condenser capacity per 24 hours	kg	2.5
Total volume of Condenser	L	4
Weight net	kg	55
Cabinet dimensions (DxWxH)	mm	500x400x520
Cabinet material		Polyester coated steel
Condenser material		Stainless steel AISI 316
Cooling media		R 507/R1150
Vacuum Valve		Solenoid
Drain tap		Yes
Keypad for control options		No
R 485/USB port for Labogene software		No

LaboGene ApS
 Bjarkesvej 5
 DK-3450 Allerød

Tel (+45) 3940 2566
 Fax (+45) 3995 2566
 Mail info@labogene.com
 Web www.labogene.com

