

characteristics Labsit 3 with glides and footring

Article: 9121 artificial leather

Dimensions in mm	Backrest:	H:	420
	Seat:	H:	520-770
		B:	450
		T:	420
Packing data	Weight	kg	9,8
	Volume	m ³	0,19



Base equipment

Benefit

Laboratory compliance	Washable und resistant to disinfectant. GMP and cleanroom compatible. Can be used in laboratories with safety classification. S1 to S3 Minimal joints and positive-locking (Hygienic-Design).	Suitable for laboratories and all hygiene-critical environments particularly health care, universities, biotechnology, chemistry and pharmaceuticals.
Handling	Practical handle on backrest. Extremely lightweight. Compact size.	Ensures maximal mobility and ease of use.
Ergonomics	Easy adjustment for user or workplace situation with single lever.	Intuitive handling. With its flexing functions Labsit adjusts automatically. No need for complicated tilt adjustments.
Backrest	Synthetic Backrest with 3D-Flex function.	The backrest adapts to the movements of the user and therefore follows the upper body when leaning back or turning sideways. This ensures a dynamic seating position and provides optimal support for the body in any position. Intuitive and automatic.
Seat	Flex funktion at front edge of seat.	The front edge of the seat automatically follows the forward-angled working position of the body. This ensures optimal blood circulation and effectively prevents pressure pain.
Seat Shell	Five colours to choose from: blue (3277), grey(3278), orange (3279) and green (3280).	Attractive design. Improves the ambience of the workplace and employee satisfaction.
Upholstered Seat	Magic Artificial Leather Washable, low maintenance, resistant to disinfectants, durable, soft and comfortable, PVC-free, phthalate-free, resistant to blood and urine, fire protection according to EN 1021 Part 1 + 2 Colour: black.	The soft, comfortable upholstery provides good seating comfort. Washable and disinfectable. Free from plasticizers so lasting softness and comfort.
Seat Height	Seat height infinitely adjustable between 520 - 770 mm.	Large variability of seat height adjustment.
5-Star Base	Sturdy polyamide base Optional: Aluminium base	Very sturdy, durable and stable. High quality. Attractive design.
Foot Ring	Height adjustable, chromium-plated foot ring.	Wide foot plate provides good support. Height can be quickly and easily adjusted.
Glides	Large, abrasion resistant floor glides.	Good stability.

characteristics Labsit 3 with glides and footing

Article: 9121 artificial leather

Dimensions in mm	Backrest:	H:	420
	Seat:	H:	520-770
		B:	450
		T:	420
Packing data	Weight	kg	9,8
	Volume	m ³	0,19



Standards	Suitable for cleanrooms category 3 in accordance with DIN EN ISO 14644-1 GMP-Conform Laboratory safety classification S1 to S3 Tested for harmful substances technical standard DIN 68 877-1 GS-Certified for highest levels of safety Blue Angel environmental certificate Sanitized Qualitymanagement DIN EN ISO 9001 Environmental Management DIN EN ISO 14001 Verified Environmental Management EMAS Occupational Health & Safety OHSAS 18001 EcoVadis CSR Rating 2016 Gold
Guarantee	10-year long-term-guarantee - 3 years of which full guarantee.

Certified by independent testing institutes. Procurement security.

Quality assured and cost effective.

characteristics Labsit 3 with glides and footring

Article: 9121 artificial leather

Dimensions in mm	Backrest:	H:	420
	Seat:	H:	520-770
		B:	450
		T:	420
Packing data	Weight	kg	9,8
	Volume	m ³	0,19



Options	Benefit
Accessories/Lap-clip	For simple and individual labelling. Easy to add or remove at any time.
Armrests	Fixed armrests. Plastic.