



## S 50 N - G 45 F Dispersing tool

/// Data Sheet

Suitable for dispersing instrument T 50 digital

Working range: 0,25 - 10 l

Stator diameter: 45 mm

Immersion depth: 70 - 250 mm

Materials in contact with medium: PTFE, AISI 316L

Suitable for solvents: yes

Max. temperature: 180 °C

Sterilization methods: all methods

## Technical Data

Volume range (H <sub>2</sub> O) [l]	0.25 - 10
Stator diameter [mm]	45
Rotor diameter [mm]	40
Gap between rotor and stator [mm]	0.5
Circumferential speed max. [m/s]	20.9
Immersion depth [mm]	70 - 250
Shaft length [mm]	290
Material in contact with medium	PTFE, AISI 316L
pH range	2 - 13
Suitable for solvents	yes
Suitable for abrasive substances	yes
Working temperature max. [°C]	180
Sterilization methods	all methods
Ultimate fineness, suspensions [µm]	10 - 30
Ultimate fineness, emulsions [µm]	1 - 10
Weight [kg]	0.929
Permissible ambient temperature [°C]	5 - 40