

## Sonifier Cell Disruptor

***Support Stand with Rod***

The support stand with stainless steel rod accommodates Branson's converter and accessories. The converter clamp supports the ultrasonic stack (converter and horn) and is easily adjusted to properly position the horn in the sample.

The stand features a black enameled cast-iron base with zinc-plated rod. A full acoustic sound enclosure is also available for sound-sensitive environments.

**STAND MEASUREMENTS:**

---

**Stand Base:** 9" l x 5.5"w x 24" h (229 x 140 mm)

---

**Rod:** 24"h x 1/2" dia. (609mm x 13mm)

---

**Part No. 101-063-1108:** Stand with rod

---

**Part No. 101-063-1109:** Converter clamp for use with Branson SFX150, SLP Sonifiers

---

**Part No. 101-063-1110:** Converter clamp for use with Branson SFX250, SFX550, 250A, 250D, 450A, 450D Sonifiers

---

# Sonifier Cell Disruptor

## Global Technologies. Local Solutions.

Branson's unmatched global resources ensure optimal solutions for the most challenging ultrasonic precision cleaning and materials joining applications. With 1600+ employees in over 70 sites worldwide, we can rapidly respond to our customers' needs, wherever they are located.

Branson understands local markets and regulations, and the open collaboration among our global offices and extensive staff of application specialists quickly produces solutions for customers. We partner with companies of

every size and scope to help resolve critical issues ranging from market changes, product quality/life cycles and production costs, to employee safety and environmental compliance. And, Branson's commitment to the publicly funded industrial collective research initiative and cooperation with various research centers have produced widespread innovations in ultrasonic cleaning and plastics joining technology that are benefiting companies everywhere.

Branson is a leading innovator in Emerson, a diversified global

manufacturing and technology company. Emerson is a progressive company working to develop products and processes to resolve global issues, including energy supply and distribution, information and networking, climate and environment, manufacturing efficiency, and product performance. Emerson is dedicated to the pursuit of breakthrough technology developments never before envisioned.

Branson brings exceptional quality, fast delivery and competitive prices. Contact your nearest regional center today.

### Americas

Branson Ultrasonics Corp.  
41 Eagle Road  
Danbury, CT 06813-1961  
T: 800-796-0400  
F: 203-796-9838  
[www.bransonultrasonics.com](http://www.bransonultrasonics.com)

### Europe

Branson Ultraschall  
Niederlassung der Emerson  
Technologies GmbH & Co. OHG  
Waldstrasse 53-55  
63128 Dietzenbach, Germany  
T: +49-6074-497-0  
F: +49-6074-497-199  
[www.branson.eu](http://www.branson.eu)

### Asia

Branson Ultrasonics (Shanghai) Co., Ltd.  
758 Rong Le Dong Road  
Song Jiang, Shanghai, PRC, 201613  
T: 86-21-3781-0588  
F: 86-21-5774-5100  
[www.branson.com.cn](http://www.branson.com.cn)