

## INC-200 Series Clear Incubators

- Clear acrylic construction allows full visibility of samples from all sides
- Choice of analog or digital temperature control
- Easy access to working chamber
- Suitable for multiple applications
- 3-year warranty



INC-200 Analog Incubator



INC-200 Digital Incubator  
shown with Rack System 01110-98

## INC-200 Series Clear Incubators

The INC-200 incubators are suitable for multiple applications such as plant propagation, humidity testing, climactic condition simulation, corrosion testing and more. The clear acrylic top and sides make it easy to view your sample from all sides.

Forced air circulation and electronic control ensures accurate temperature conditions. A pre-set push switch automatically sets the temperature at 37 °C, while a safety feature cuts off heating if the unit reaches 72 °C or higher.

Temperature is set via a calibrated knob for the analog model (INC-200), and via an LCD when using the digital model (INC-200D). The latter not only allows more accurate temperature management but also gives a constant reading of actual temperature when in use.

The purpose-built design offers easy access to samples – a hinged front door panel and plugged holes (2 x 10 mm diameter) on each of the side panels permits the introduction of gases or cables. These units do not have a base so they can be placed directly over entire instruments; order optional base if needed.

### Specifications

Model no.	INC-200
Temperature range	Ambient +5 °C to 60 °C
Temperature fluctuation at 37 °C	±0.1 °C
Temperature variation between shelves	±0.3 °C
Capacity	60 L (2.1 cu ft)
Internal dimensions (W x H x D)	45 x 38 x 38 cm (17.7 x 15 x 15 in)
Overall dimensions (W x H x D)	60 x 39 x 39 cm (23.6 x 15.4 x 15.4 in)
IP rating	IP30
Warranty	3 years

### Ordering Information

Description	Item no.	Series no.	Model no.	Legacy SKU
Clear analog incubator, 230 VAC, 50 Hz	01110-02	INC-200	INC-200	SI60
Clear analog incubator, 120 VAC, 60 Hz	01110-00	INC-200	INC-200-120	SI60/120V/60
Clear digital incubator, 230 VAC, 50 Hz	01110-08	INC-200	INC-200D	SI60D
Clear digital incubator, 120 VAC, 60 Hz	01110-04	INC-200	INC-200D-120	SI60D/120V/60

### Accessories

Description	Item no.	Series no.	Model no.	Legacy SKU
Acrylic base plate	01110-97	INC-200	N/A	5160/1
Plastic coated rack system with two shelves	01110-98	INC-200	N/A	SI60/2
Sample tray	01110-99	INC-200	N/A	SI60/5

coleparmer.com

USA:  
+1.800.323.4340  
+1.847.549.7600

Canada: +1 800 363 5900  
China: +86 21 5109 9909  
France: +33 (0)1 48 63 78 00  
Germany: +49 (0) 9377 92030

India: +1 800 266 1244  
Italy: +39 (0)2 84349215  
UK: +44 (0) 1480 272279  
All other countries: +1 847 549 7600

**Cole-Parmer®**