

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Laminin from Engelbreth-Holm-Swarm murine sarcoma basement membrane - 1-2 mg/mL in Tris buffered NaCl, 0.2 µm filtered, BioReagent, suitable for cell culture

Product Number:

L2020

CAS Number:

114956-81-9

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Form)	Solution
Appearance (Colour)	Colorless
Appearance (Turbidity)	Clear to Slightly Hazy
Attachment (Attach./Toxicity)	Pass
Protein Content (Method)	Lowry
Prot. Content (mg Protein/ml)	1.0 - 2.0
Vial Content (Volume ml)	0.5 - 1.0 ml
Vial Content (Size)	1.0 - 1.0
Vial Content (Unit Of Measure)	mg
Vial Content (Actual Content)	1.0 - 1.3
Vial Content (Overpack)	≤ 30.00 %
Microbial Contamination	Pass
Gel Electroph. (Gel)	SDS-Page
Gel Electroph. (Stain)	Coomassie
Gel Elec. (Compar. to Cont)	Similar
Expiration Date Period	-----
2 Years	

Specification: PRE.3.ZQ5.10000003248

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.