

# Cole-Parmer®

## Cole-Parmer® Slide Drying Bench - SW-200

This Slide Drying Bench is used for the safe drying of glass slides after sample mounting in preparation for clearing and staining. The dark surface gives good visibility of samples.

- Up to 50 slides capacity
- The slides can be arranged on the bench in different orientations
- Built-in energy regulator controlling temperature up to 100°C



# Cole-Parmer® Slide Drying Bench - SW-200

Cole-Parmer® SW-200 drying bench aids in the preparation of microscope slides at the specimen mounting stage. It complements both the paraffin wax dispenser and the paraffin wax mounting bath to provide a complete set of paraffin wax sampling equipment for the modern laboratory.

Accepting up to 50 slides, it has the facility for drying slides in different orientations- across the drying support bars, angled from the bars, or flat without the bars.



## Technical Specification

Specification	SW-200
Capacity	Up to 50 slides
Material	Base and top powdercoated
Max Element Temperature	100°C
Dimensions (d x w x h)	180 x 390 x 95mm
Weight	1.8kg

## Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Microscope Slide Warmer/Dryer 230V UK & EU	48573-06	SW-200	SW-200	MH6616
Microscope Slide Warmer/Dryer 115V US	48573-01	SW-200	SW-200-115	MH6616X1

[coleparmer.com](http://coleparmer.com)

USA:  
+1.800.323.4340  
+1.847.549.7600

Canada: +1.800.363.5900  
China: +86.21.5109.9909  
France: +33 (0) 1486 37800  
Germany: +49 (0) 9377 92030

India: +1.800.266.1244  
Italy: +39 (0) 1313 89513  
UK: +44 (0) 1480 272279  
All other countries: +1.847.549.7600

**Cole-Parmer®**  
essentials  
an antylia scientific company