

IKA

designed for scientists



MultiDrive BL 2000.1 Standard blend knife, stainless steel

/// Data Sheet

Standard mixing knife, stainless steel

Compatible with IKA MultiDrive BL 2000 / 2000 T



designed for scientists

Technical Data

Process type	batch
Operating principle	cutting
Feed hardness max. [Mohs]	5
Material beater/cutter	stainless steel 1.4034
FDA conform	yes
Multi use	yes
Beater/Knife diameter [mm]	75
Permissible ambient temperature [°C]	5 - 40