



designed for scientists



## MultiDrive MI 400.1 Standard beater, stainless steel

/// Data Sheet

Standard beater, stainless steel

Suitable for most brittle materials up to Mohs hardness 5

Compatible with IKA MultiDrive MI 400 / 400 T



designed for scientists

### Technical Data

Process type	batch
Operating principle	impact
Feed hardness max. [Mohs]	5
Material beater/cutter	stainless steel 1.4034
FDA conform	yes
Multi use	yes
Beater/Knife diameter [mm]	81
Permissible ambient temperature [°C]	5 - 40