



designed for scientists



MultiDrive MT 150.5 Milling tube, autoclavable

/// Data Sheet

MT 150.5 is a plastic milling tube with a stainless steel beater, which is ideally suited for grinding dry samples in combination with the MultiDrive control universal crusher. Thanks to the transparent plastic, the grinding material can be continuously observed.

An integrated RFID chip simultaneously stores and limits the lifespan of the milling tube. 5 spare seals are included in the scope of delivery. By replacing the seal, the service life of the grinding chamber can be extended. Cross contamination is effectively prevented by sterilizing the milling tube in an autoclave.



designed for scientists

Quantity: 5 tubes/ pack

Useable volume max.: 150 ml

Feed grain size max.: 10 mm/ 5 Mohs

Material in contact with medium: Stainless steel, PA, PC, PTFE, silicone

FDA compliance: Yes

Note: For safety reasons TC 1 protective housing is mandatory to operate the tube. MultiDrive won't start, if TC 1 isn't attached. TC 1 is available as an accessory.



IKA

designed for scientists

