



designed for scientists



## MultiDrive BT 250.5 Blending tube, autoclavable

/// Data Sheet

BT 250.5 is a plastic blending tube with a stainless steel beater, which is ideally suited for wet grinding in combination with the MultiDrive control universal crusher. Thanks to the transparent plastic, the change of the sample material can be continuously observed.

An integrated RFID chip simultaneously stores and limits the lifespan of the blending tube. 5 spare seals are included in the scope of delivery. By replacing the seal, the service life of the grinding chamber can be extended. Cross contamination is effectively prevented by sterilizing the grinding chamber in an autoclave.



designed for scientists

Quantity: 5 tubes/ pack

Useable volume max.: 250 ml

Feed grain size max.: 10 mm/ 5 Mohs

Material in contact with medium: Stainless steel, PA, PC, PTFE, silicone

FDA compliance: Yes

Note: For safety reasons TC 1 protective housing is mandatory to operate the tube. MultiDrive won't start, if TC 1 isn't attached. TC 1 is available as an accessory.





designed for scientists

