

Cole-Parmer®

Cole-Parmer® See-saw Rockers - SH-200D-S, SH-200D-S-L

- Small space saving design - ideal for personal use
- Gentle See-saw rocking action
- Incubator safe
- Digital speed control and built-in timer
- Optional tier system available to increase capacity



SH-200D-S
Mini See-saw Rocker



SH-200D-S-L See-saw Rocker

Cole-Parmer® See-saw Rockers - SH-200D-S, SH-200D-S-L

- Small space saving design or laboratory sized unit
- Digital speed control and built-in timer
- Optional tier system available to increase capacity

The Cole-Parmer® models have a see-saw rocking action that creates a wave motion within vessels such as culture flasks, Petri dishes etc. Digital selection of both speed and time facilitates accurate and reproducible conditions. Both parameters are controlled via an easy to use control dial and displayed on the bright LED display.

Both models are supplied with a non-slip mat which holds flat based vessels securely in place during mixing. An accessory tier system can be fitted in seconds, without the need for tools, tripling the available space for samples. Each tier is securely held in place by magnets.

Technical Specification

Specification	SH-200D-S	SH-200D-S-L
Rocking action	See-saw wave	See-saw wave
Platform dimensions (w x l)	235 x 235mm	355 x 355mm
Speed range	5 to 70rpm	5 to 70rpm
Angle of tilt	7°	7°
Maximum load	3kg	10kg
Overall dimensions (w x d x h)	240 x 300 x 150	360 x 420 x 170
Tier height, mm	125/tier	125/tier
Operational temp. range	+4 to +40°C	+4 to +40°C
Max. permissible humidity	80%	80%
Net weight	5kg	10kg
Electrical supply	120V, 60 Hz, 50W 230V, 50Hz, 50W	120V, 60 Hz, 50W 230V, 50Hz, 50W

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Cole-Parmer® See-saw rocker mini	51901-14	SH-200	SH-200D-S	SSM4
Cole-Parmer® See-saw rocker mini 120V	51901-15	SH-200	SH-200D-S-120	SSM4/120V/60
Cole-Parmer® See-saw rocker large	51901-06	SH-200	SH-200D-S-L	SSL4
Cole-Parmer® See-saw rocker large 120V	51901-07	SH-200	SH-200D-S-L-120	SSL4/120V/60

Accessories

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Cole-Parmer® Tier system (2 platforms plus 8 bars)	51700-93	SH-200	SH-200D-S/SH200D-S/120	SSM3/1
Cole-Parmer® Tier system (2 platforms plus 8 bars)	51700-94	SH-200	SH-200D-S-L/SH200D-S-L/120	SSL3/1

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

Cole-Parmer®
essentials

an antylia scientific company

Find out more at antylia.com