


This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Sign reveal of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.



ISSUE	DATE	AMENDMENT	ISSUE	DATE	AMENDMENT



 DO NOT SCALE IF IN DOUBT ASK	CUSTOMER CUSTOMER DRG No MATERIAL FRAME: Polycarbonate, TUBES: Polypropylene TOLERANCES GENERAL TOLERANCES (UNLESS OTHERWISE STATED) <small>0 ±0.20 0.05 ±0.10 0.10 ±0.15 ANGULAR ±0.30 (ALL DIMENSIONS IN MILLIMETRES)</small>	DRAWN ZOLTAN E DATE 31.10.2014 CHECKED JAMES B DATE SCALE	TITLE Vari-Strip™ clear, Pre-loaded in fully skirted frame Marketing drawing DRAWING No 4ti-0753/C/757 SHEET No 1 OF 1 ISSUE A
-------------------------------------	--	---	--