

**Product Name: UniRack™**

**Catalogue No.: S500-80 series**

**This document replaces any previous version**

## 1. Product Description:

- **Rack:** Non-sterile, 80 position rack
  
- **Available colors:**
  - Blue, Red and Yellow
  
- **Compatible with:**
  - 1 to 5 ml Cryotubes and 0.5 to 2 ml microtubes
  - 0.2 to 0.5 ml PCR tubes

## 2. Packaging:

- **Case:** 10 units per case

## 3. Product Specifications:

- **Material:**
  - Polypropylene
- **Temperature range:** -196°C to +121°C
- **Autoclavable** at 121°C, for up to 30 minutes
- **Total capacity:** 80 cryovials and microtubes on one side or 60 PCR tubes on other side

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;

## 8. Recommended Use:

Rack for holding tubes

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*