


Fully assembled, mega racks are great for high throughput labs requiring large sample storage

- Sturdy, durable material for easy and secure transport of clinical and research tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Reusable, easy to clean
- Place different test groupings into one rack eliminating the need for multiple racks
- Stackable for staging and saving space
- Ideal footprint for use with programmable tube filling machines




120114

Mega Racks® 10-13 mm - Holds 5-7 mL sample tubes



Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Material
HS24322A	Blue ■	216	13	13.3 x 9.3 x 3.6	33.7 x 23.5 x 9.1	PP
HS24322E	White □					
120112	White □	432	13	26.5 x 9.3 x 3.5	67.3 x 23.5 x 9.0	PP
120113	Blue ■					



HS24323A & HS24323E

Mega Racks® 13-16 mm - Holds 10 mL sample tubes



Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Material
HS24323A	Blue ■	120	16	13.3 x 10.5 x 3.6	33.7 x 26.7 x 9.1	PP
HS24323E	White □					
120114	White □	240	16	26.5 x 10.5 x 3.5	67.3 x 26.7 x 9.0	PP
120115	Blue ■					

Mega Racks® 18-20 mm - Holds 15 mL conical/swab tubes and 20 mm round bottom tubes



Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Material
120834	White □	120	20	14.9 x 11.9 x 3.6	37.7 x 30.1 x 9.1	PP
120835	Blue ■					




120834

Mega Racks® 28-32 mm - Holds 50 mL conical/swab tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Material
120958	Blue ■	100	32	16.2 x 16.7 x 4.7	41.2 x 42.3 x 12	ABS
*120959	White □					

*Not recommended for autoclaving



120959