

Product Name CELL156



Key Features

- +2°C to +10°C temperature range
- Annual service & Calibration available
- Fitted lock with 2 keys



Flexible Storage

Multiple shelf locations. Back stop on every shelf provides a flexible storage solution.



Spark Free

No internal light or fan, suitable for a spark free laboratory.

Product Name CELL156

Details	
Net Capacity (litres)	156
Unit External Width (mm)	600
Unit External Depth (mm)	600
Unit External Height (mm)	835
Product Weight (kg)	34
Packed Width (mm)	630
Packed Depth (mm)	644
Packed Height (mm)	885
Packed Weight (Kg)	36
Temperature Range	+2°C to +10°C
External Finish	Painted Steel
Internal Finish	White Plastic
Battery Back up	No
Door Finish	Painted Steel
Door Type	Solid
Door Style	Hinged, Right hand - Reversible
Door Open Depth (mm)	1150
Replacement Door Seal	No
Open Door Sensor	No
Lock	Yes
Adjustable Feet (number)	Yes, 2
Porthole	No
Porthole location	N/A
Temperature Display on Panel	No
Temperature Alarm	No
Remote Alarm Contacts	No
Audible Alarms	No
Visual Alarms	No
Number of Shelves	3
Number of Shelf Positions	5
Shelf Material	Glass
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	No
Spark Free	Yes
Controller Type	Manual
Number of Probes	1
Probe Type	Air
Ambient Temperature Range	+10c to +30c
Refrigerant Type	R600a
CFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Voltage Supply	220-240V / 50Hz
Certification	CE/UK
Noise Level (dba)	41
Energy Consumption (KWh/annum)	117
Energy Rating	F