

# LGPv 1420

## MediLine



<b>Total gross capacity</b>	1427 litre
<b>Ext. dimensions</b> in mm (w/d/h)	1430/830/2150
<b>Int. dimensions</b> in mm (w/d/h)	1250/700/1550
<b>Energy efficiency class</b>	
Energieverbrauch in 365 Tagen	3244 kWh
Climate rating <sup>1</sup>	5
Cooling system	–
Defrosting	automatic
Temperature range	-10°C to -26°C
Gastronorm	- mm
Housing material/colour	
Door/lid material	Steel
Liner Material	304-grade stainless steel
Type of control	electronic with text display
Temperature display	external digital
Interior light	–
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Shelves	8
Adjustable Shelving	–
Number of drawers	–
Number of baskets	–
Maximum number of baskets	–
Shelf material	grid shelves plastic-coated
Shelf loading	60 kg
Height of adjustable legs	-
Castors/type of castors	yes/swivel castors
Door handle	full-length integrated door handle
Foot pedal opener	–
Self-closing door	yes
Insulation	83 - 60 mm
Defrost water drain	–
Refrigerant	R 290
Airborne acoustical noise emissions	60 dB(A)
Gross/net weight	262,0/219,0 kg
EAN-Nr.	9005382137215

1 Climate Ratings: **SN** +10°C to +32°C    **N** +16°C to +32°C    **ST** +16°C to +38°C    **T** +16°C to +43°C