



**STANDARD series**

## Upright Refrigerator Model: BIO360FRS

The **BioCOLD STANDARD** series of refrigeration has been built to a very high specification, to meet the requirement for general laboratory storage.

The Upright Refrigerator features 5 glass safety shelves as standard.



### Standard Features

- Glass shelves
- Reversible doors
- Replaceable door gasket
- CFC and HCFC free refrigeration system

### Options

- Sparkfree interior
- Door mounted temperature display/alarm module

### Technical Specification

#### Dimensions

External (H x W x D), (mm)  
Weight (kg)

1750 x 600 x 655  
75

#### Temperature

Set Point (adjustable) (° C)  
Operating range (° C)

+4  
+2 to +8

#### Capacity

Refrigerator Volume (L)

350

#### Electrical Supply

Voltage (V)  
Frequency (Hz)  
Fuse (A)

220-240  
50  
13

#### Shelves (Glass)

5

#### Defrost Type

Manual

**BioCOLD**

Telephone/Fax: 01506 873795 E-mail: [sales@biocold.co.uk](mailto:sales@biocold.co.uk) [www.biocold.co.uk](http://www.biocold.co.uk)