



**STANDARD series**

## Upright -20°C Freezer Model: BIO270FZSS

The **BioCOLD STANDARD** series of refrigeration has been built to a very high specification, to meet the requirement for general laboratory storage.

The Upright -20°C Freezer features 2 drawer flaps and 5 removable solid baskets.



### Standard Features

- Sparkfree interior
- High temperature alarm (visual)
- Reversible doors
- Replaceable door gasket
- CFC and HCFC free refrigeration system

### Options

- Door mounted temperature display/alarm module
- Door lock

### Technical Specification

**Dimensions**

External (H x W x D), (mm)	1750 x 600 x 600
Weight (kg)	63

**Capacity**

Freezer Gross, (L)	270
--------------------	-----

**Drawers** (solid baskets)

2/5
-----

**Defrost Type**

Manual
--------

**Temperature**

Set Point (adjustable) (°C)	-20
Operating range (°C)	-15 to -25

**Electrical Supply**

Voltage (V)	220-240
Frequency (Hz)	50
Fuse (A)	13



**Telephone/Fax:** 01506 873795 **E-mail:** [sales@biocold.co.uk](mailto:sales@biocold.co.uk) **www.biocold.co.uk**