


FREE STANDING

# Solid Door: DLS700PS

## Overview

Our intuitive range refrigerators feature built-in wireless monitoring software, enabling you to access the refrigerator's performance remotely.

Plus, our bespoke touchscreen controller, designed specifically for the end user with simple access to key information.



## Key Features



### Advanced Controller

Easy to use touch screen controller, with no awkward buttons to press



### Connected

Wireless temperature monitoring built-in



### Dual Temperature Probes

View your daily minimum and maximum temperatures in one touch



### Battery Back-Up

A/V alarms should your temperature stray or the door be left open



### Antimicrobial Protection

Biocote® handle helping protect from microbes such as E. coli, Campylobacter, Salmonella and MRSA



### Self-Closing Doors

Fitted with self-closing hinges to help avoid any temperature loss

# DLS700PS - Product Specification

| Counter Top                      |                          |
|----------------------------------|--------------------------|
| Model                            | DLS700PS                 |
| Gross Capacity (Litres)          | 700                      |
| Net Capacity (Litres)            | 482                      |
| Unit External Width (mm)         | 720                      |
| Unit External Depth (mm)         | 800 (excl handle)        |
| Unit External Height (mm)        | 2040                     |
| Packed Width (mm)                | 810                      |
| Packed Depth (mm)                | 910                      |
| Packed Height (mm)               | 2210                     |
| Packed Weight (kg)               | 125                      |
| Internal Dimensions (mm) (WxDxH) | 600x665x1530             |
| Temperature Range                | +2C to +8C               |
| External Finish                  | Pre-coated steel         |
| Internal Finish                  | Stainless Steel          |
| Type of Plastic                  | N/A                      |
| Door Type                        | Solid / Self Closing     |
| Door Style                       | Reversible               |
| Open Door Depth                  | 1465mm                   |
| Replaceable Door Seal            | Yes                      |
| Open Door Sensor                 | Yes                      |
| Lock                             | Yes / Mechanical         |
| Adjustable Feet                  | Castors x4               |
| Porthole                         | Yes                      |
| Porthole Size                    | 32mm                     |
| Temperature Display On Panel     | LED / Touchscreen        |
| Temperature Alarm                | Yes                      |
| Remote Alarm Contacts            | Yes                      |
| Audible Alarms                   | Yes                      |
| Visual Alarms                    | Yes                      |
| Number Of Shelves                | 5                        |
| Number Of Shelf Positions        | 20                       |
| Shelf Size                       | 530x650mm                |
| Shelf Material                   | PVC coated wire          |
| Number Of Baskets                | 0                        |
| Basket Size                      | N/A                      |
| Basket Material                  | N/A                      |
| Interior Light                   | Yes                      |
| Light Type                       | LED                      |
| Light Activation                 | Door Switch / Controller |
| Fan                              | Yes                      |
| Spark Free                       | No                       |
| Defrost Setting                  | Auto                     |
| Controller Type                  | Digital / Wireless       |
| Battery Back Up                  | Yes                      |
| Controller Connectivity          | Wireless / SD Download   |
| Number Of Probes                 | 2                        |
| Probe Type                       | Air / Load               |
| Ambient Temperature Range        | 10 - 32C                 |
| Refrigerant Type                 | R290                     |
| Plug Type                        | UK                       |
| Voltage Supply                   | 220-240V                 |
| Frequency (Hz)                   | 50Hz                     |
| Insulation Thickness             | 58mm                     |
| Anti Microbial Handle            | No                       |