



NINKA

STYLE & PERFORMANCE

Protective eye shield cost effective and accessible to all that was specifically designed for the service and healthcare industries. The medical and retail sectors will also appreciate its ergonomics and protection against splashes. Light yet resistant, NINKA ensures a safer work environment for every person who may be in close physical contact with others.

- + **ERGONOMIC TIPGRIP TEMPLES WITHOUT PRESSURE POINTS**
- + **SIMPLE & QUICK LENS REPLACEMENT**
- + **PROTECTION AGAINST SPLASHES & DROPLETS**
- + **UPPER FRAME PROTECTION**



TECHNOLOGIES

VERSIONS



PSONINK010



NINKA

VERSION	REFERENCE	LENS MARKING	FRAME MARKING	COATING	STANDARD
CLEAR	PERSONINK010	CE	CE	ANTI-FOG	CE

BOLLÉ SAFETY PRODUCT APPROVAL CERTIFICATES TO CE STANDARDS WILL BE SENT ON REQUEST.


PACKAGING & EAN13 CODE

SINGLE:	INNER BOX:	OUTER PACK:
Size: 38 x 14.5 x 19 mm	Size: xx cm Weight: g	Size: 60 x 40 x 30 cm Weight: 10.7 kg
3660740010690	- pcs	3 660 740 200 695 -20 pcs
CUSTOM CODE: 9004901000		

GENERAL SPECIFICATIONS

LENS	PET	STRAP	
FRAME	PP - HALF FRAME	USE	SINGLE USE
WEIGHT	1610 g	SUPPLIED WITH	

KNOW MORE
ABOUT
MARKINGS:



OR SCAN:

CLEANING & STORAGE

CLEANING: The lenses should be cleaned regularly. Bollé Safety recommends only using the B-clean lens cleaning system.

MAINTENANCE: Pitted, scratched or damaged lenses must be replaced by replacing the complete product. If the frame is damaged also replace the complete product. After each use it is recommended to keep the spectacles or goggles in a microfiber pouch or rigid case.

STORAGE: All products must be stored in a dry and clean environment.



MADE IN: TW

DATASHEET: 01-07-2022

bolle-safety.com

Europe Headquarters
34 rue de la Soie
69100 Villeurbanne
France
Tél. : 00 33 (0) 4 78 85 23 64

U.K. & Ireland
Unit C83 - Barwell Business Park
Leatherhead Road - Chessington
Surrey - KT9 2NY - UK
Tél. : 00 44 (0)208 391 4700

USA Office
2810 Caribou Court, Suite 160
Carlsbad, CA 92010
United States of America
Tél. : 800 222 6553

Asia-Pacific Office
32-40 Fairchild Street
heatherton VIC 3202
Australia
Tél. : 1300 138 098



**THE PPE EYEWEAR
SPECIALISTS**