

# KIMTECH™

## Kimtech™ Polaris™ Xtra Nitrile Gloves

- Textured fingertips
- Enhanced chemical and biohazard protection
- Optimal skin health protection

### BUILT FOR SCIENCE – MADE FOR YOU

Kimtech™ Polaris™ Xtra Nitrile Gloves have been designed for use in laboratories, life science applications and non-sterile pharmaceutical manufacturing environments.

With extend cuff length and proven protection against the harshest chemicals, superior dexterity and comfort you can count on, these gloves are built for science, made for you.



CE 2777

UK CA 0321

- +29 common chemicals tested
- +24 Cytotoxic drugs tested

Recyclable with our RightCycle™ Programme



Size Guide	XS (5.5)	S (6.5)	M (7.5)	L (8.5)	XL (9.0)
Codes	62100	62101	62102	62103	62104
Packaging	50 gloves per dispenser - 10 dispensers per case (=500 gloves)				

## Description

<b>Formulation</b>	High grade synthetic nitrile rubber ( acrylonitrile - butadiene)
<b>Technology</b>	Optimized nitrile formulation for extended chemical splash protection with no compromise on wearer comfort
<b>Design</b>	Dark Magenta, Ambidextrous, beaded cuff, textured fingertips

## Assured Compliance

<b>CE/UKCA</b>	PPE Cat III according to Regulation (EU) 2016/425 and to the Regulation 2016/425 as brought into UK law and amended.
<b>EU PPE norms</b>	EN ISO 374-1 Type B (JKOPT) Chemical Splash Protection EN 374-4 Resistance to Degradation by Chemicals EN 374-5 Micro-Organism and Virus Protection
<b>Other Standards</b>	ASTM D6978 Chemotherapy Drug Splash Protection ISO 10993 Bio-compatibility Food contact approved. Compliant with German BfR Recommendation XXI and XXI/1 and Regulation (EU) 2019/515
<b>Quality Assurance</b>	Quality System of the Manufacturing Facility compliant with ISO 9001, ISO 13485 and 21 CFR Part 820

## Residual Content

<b>Vulcanization Agent</b>	Non detectable level of vulcanization agent in the finished product
<b>Powder</b>	Powder-free glove with a maximal residual powder content <2.0 mg/glove
<b>Latex Protein</b>	Not made with natural rubber latex

## Documentation

<b>Certificate of Analysis (CoA)</b>	Certificate of Analysis available for each production lot at <a href="http://www.certs.kcprofessional.com">www.certs.kcprofessional.com</a>
<b>Declaration of Conformity (DoC)</b>	Declaration of Conformity available from the product page on <a href="http://www.kimtech.eu">www.kimtech.eu</a>
<b>User Instructions</b>	Printed on the dispenser box

## Product Specifications

Characteristics	Value					Test method
Freedom from holes	AQL 0.65 (as defined per ISO 2859-1 sampling by attributes)					EN 374-2
<b>Tensile properties</b>	Tensile Strength (min.)		Elongation at Break (min.)			
Before Aging After Aging	≥ 14 MPa ≥ 14 MPa		≥ 500%, min. ≥ 400%, min.			ASTM D412 ASTM D573 ASTM D6319
<b>Dimension</b>	<b>Target Thickness/Width/ Length</b>					
Thickness (mm)	Middle Finger 0.19	Palm 0.14		Cuff 0.11		ASTM D6319 EN 455-2
Palm width (mm)	XS 70	S 80	M 95	L 110	XL 120	ASTM D6319 EN 455-2 EN ISO 21420
Length (mm)	300					

Visit us at [www.kimtech.eu](http://www.kimtech.eu) or for any questions, email: [kimtech.support@kcc.com](mailto:kimtech.support@kcc.com)



RightCycle by Kimberly-Clark Professional® is an innovative recycling programme that diverts non-hazardous Kimtech™ branded nitrile gloves and coveralls from landfill. The programme is available in the UK, Ireland, France, Spain, Germany, Austria, Belgium, Netherlands and Switzerland.

The color Dark Magenta is a trademark of Kimberly-Clark Worldwide Inc.

©™ Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. © KCWW. Publication code: ID 5466 EN 01.24

