

# DFS2

**DFS2** is a **facial screen** accessible to all which perfectly protects against splashes and droplets. **Designed to fit the first responders' needs**, DFS2 is simple to use, economical, protective and can be worn over glasses.




PFSDFS2106



52 g

- + **READY-TO-USE PACKAGING**
- + **ANTI-FOG MATERIAL & PROTECTIVE FILM ON OUTER SIDE**
- + **ADAPTED TO PRESCRIPTION GLASSES**
- + **PROTECTION AGAINST SPLASHES & FLUIDS**
- + **ELASTIC HEADBAND (ONE SIZE FIT ALL)**
- + **XL SIZE VISOR: 22 X 33CM**
- + **EN166: PPE CATEGORY II**

# DFS2

MODEL	VERSIONS	REFERENCES	LENS MARKINGS	FRAME MARKINGS	COATINGS
DFS2	 Clear	PFSDFS2106	CE	-	AF

## CLEANING AND STORAGE

**CLEANING:** The lenses should be cleaned regularly.

**MAINTENANCE:** Pitted, scratched or damaged lenses must be replaced. If the frame is damaged also replace the product.

**STORAGE:** All products must be stored in a dry and clean environment.

**PACKAGING:** Each single packaging is sealed and contains user instruction.

## LENS MARKINGS

CE European standard

Made in China

## PACKAGING AND EAN13 CODE

	PFSDFS2106
Plastic bag of 5 products	3660740010577
Outer Pack of 200 products	3660740200572
Outer pack size: 69 x 38 x 35cm Outer pack G.W : 12.15 Kg	

## CHARACTERISTICS

VISOR	PET - 0.3 mm
FRAME	PU (no latex)
STRAP	White Polyester
SPARE LENS	No
PROTECTION	PPE CAT 2 COVID

01/04/2020

EMEA  
Bollé Safety  
34 rue de la Soie  
69100 Villeurbanne  
France

[bolle-safety.com](http://bolle-safety.com)

**bolle**  
SAFETY