



Cleanroom Solutions Engineered for Protection, Compliance and Comfort

The Kimtech™ A5 Sterile Cleanroom Coveralls simplify aseptic gowning to save time and reduce the risk of contamination



CLEAN-DON™ blue line on the inside indicates correct place to grasp for aseptic gowning.



Telescope fold and snaps prevent legs and arms falling to the floor during gowning.



Vacuum packed and double bagged for added sterility assurance. Reduces space so more products fit into gowning area storage.

The Kimtech™ A5 Sterile Integrated Hood & Mask is designed to guard the face and head from common contaminants, providing seamless protection



No gaps between hood and mask ensures exhaled air is filtered. **Elastics** provide a comfortable, secure and universal fit.



CLEAN-DON™ blue line on the inside indicates correct place to grasp for aseptic gowning. Saves time and reduces contamination risk.



Vacuum packed for added sterility assurance. Reduces space so more products fit into gowning area storage.

The Kimtech™ A5 Sterile Boot Covers improves worker safety and comfort, while reducing the risk of contamination



New version ideal for wet applications. Wrap-around vinyl foot is seamed above the floor line to reduce risk of liquid penetration and trip hazards. Standard vinyl sole is also available.



3 sizes for better safety and comfort. Better fit and comfort to improve worker safety by reducing risk of trip hazards.



Vacuum packed for added sterility assurance. Reduces space so more products fit into gowning area storage.

Applications

EU GMP ISO 5 Grade A	Aseptic Processing	Parental Drug Manufacturing	Biotechnology	Pharmaceutical Compounding	Ophthalmic Product Manufacturing	Aseptic Hood work
-------------------------	-----------------------	--------------------------------	---------------	-------------------------------	--	----------------------

Fabric Tests	Test Method	Result
Particle shedding test (Helmke Drum)	IEST-R-CC003.3	Category I
Bacterial filtration efficiency (3.0µm)	ASTM F2100	96%
Particle filtration efficiency (3.0µm)	ASTM F2299	94%
Abrasion resistance	EN 530 Method 2	Class 1
Flex cracking resistance	ISO 7854 Method B	Class 6
Trapezoidal tear resistance	ISO 9073-4	Class 1
Puncture resistance	EN 863	Class 1
Tensile strength	EN ISO 13934-1	Class 1
Resistance to ignition	EN 13274-4 Method 3	Pass
Seam strength	EN ISO 13935-2	Class 3

Tests performed on fabric used in all Kimtech™ A5 Sterile Cleanroom Apparels.

Kimtech™ A5 Sterile Cleanroom Coveralls

Codes and Sizes	Case Count
88800 (S), 88801 (M), 88802 (L), 88803 (XL) 88804 (2XL), 88805 (3XL), 88806 (4XL)	25 coveralls per case

Kimtech™ A5 Sterile Integrated Hood & Mask

Code	Case Count
Integrated Hood & Mask (36072)	75 individual units per case
Standard Hood with Ties (25797)	100 individual units per case

Kimtech™ A5 Sterile Cleanroom Boot Covers

Code	Case Count
31683 (S), 31696 (Universal)	100 individually packed pairs per case
12922 (S), 88808 (Universal)	100 individually packed pairs per case



Visit www.kimtech.eu for more product information: includes data-packs, gowning posters and certificates of conformance & irradiation.