

Standard


Vinyl Examination Gloves



Glove Plus Vinyl examination gloves offer cost effective barrier protection.

Manufactured in a dedicated Vinyl factory, this product is 100% Natural Rubber Latex free and is ideal for glove users who have been identified as having allergic reactions to NR Latex proteins (Type 1 reactions) in areas of low risk.

Product Key Points



Product Information

Type: Powder free & non-sterile.

Material: Vinyl - Containing DINP Plasticizer

Protein Content: 100% NR Latex free.

Powder Content: Below 2mg/g glove.

Features: Ambidextrous, Beaded Cuff.

Colour: Opaque.

Packing: 100 gloves per dispenser, 10 dispensers per carton.

Dimensions

Glove Size	Product Code	Palm Width (mm)	Length (mm)
Extra Small	SV31	< 80	245
Small	SV32	± 80	245
Medium	SV33	± 95	245
Large	SV34	± 110	245
Extra Large	SV35	> 110	245

Thickness

Location of Thickness Measurements	Single Wall (mm)
Finger (at 15mm at the tip)	0.19 ± 0.03
Palm (at centre of palm)	0.17 ± 0.03

Physical Properties

Parameters	Before Aging	After Aging
Force at Break (N)	7.0	6.0
Watertight AQL	1.5	N/A

Pre-Shipment Quality Inspection

Criteria	Insp Level	AQL
Dimensions	S-2	4.0
Physical Properties	S-2	4.0
1000ml Water Leak Test	G-1	1.5
Visual Inspection	S-4	4.0

Product Conformance:

- EU MDD EN455 Parts 1 & 2.
- EN374 Parts 1, 2 & 3.
- EN420 Part 3.

Quality Assurance:

- This product is manufactured in a plant where the quality management system has been independently assessed as conforming to ISO9001:2008 and ISO13485:2003.



www.glove-plus.com