

**ID-Nr. 30495-041**  
**Model CS.195.105**



**Cabinet for chemicals CS-CLASSIC**  
**Model CS.195.105, body colour anthracite grey**  
**(RAL 7016), with wing doors in light grey (RAL 7035)**

**Safe and proper storage of non-flammable chemicals**  
**and toxic substances in working areas**

**Function / construction:**

- **Robust construction and longevity:** robust, high-quality, scratch-resistant outer body with powder-coated textured surface
- **No unauthorised use:** doors lockable with cylinder locking (integrated in the turning handle)
- **Easy alignment:** adjusting aids to compensate for uneven floor
- **Ventilation:** natural ventilation due to ventilation openings at the cabinet bottom, ready for connection to a technical exhaust system

Dimensions	
Width (external)	1055.00 mm
Depth (external)	520.00 mm
Height (external)	1950.00 mm
Width (internal)	996.00 mm
Depth (internal)	493.00 mm
Height (internal)	1835.00 mm

The available internal dimensions can vary depending on the interior equipment package.


Weight / Load	
Weight without interior equipment	75 kg
Distributed load	421.00 kg/m <sup>2</sup>
Maximum load	600 kg

Ventilation	
Extraction air	75 DN
Minimum exhaust air quantity 10 times	9 m <sup>3</sup> /h
Differential pressure (10 times air change)	2 Pa


Colouring / Surface finish	
Body colour code	RAL 7016
Body colour	anthracite grey
Body surface finish	textured
Door / Drawer colour code	RAL 7035
Door / Drawer colour	light grey
Door / Drawer surface finish	textured

General Information	
Number of doors	2
Depth with open doors	1026 mm
Customs tariff number	94032080

Front view



Side view



Top view

