

Nitrile Examination Gloves

Latex free, powder free, non-sterile

HAIKA® NX510 are lightweight, ambidextrous, single-use medical nitrile examination gloves.

They offer healthcare professionals high performance barrier protection as well as supreme comfort and touch sensitivity.

Key Features

Secure Control

A textured surface enhances slip resistance and allows for a secure grip.

Supreme Comfort and Sensitivity

The lightweight and elastic nature of the nitrile delivers an excellent sense of touch and a comfortable wearing experience even when worn for long periods.

Effective Barrier Protection

Helps to protect users and patients from the transmission of infections.

Beaded Cuff

Durable beaded cuff design enables the glove to be donned and removed efficiently and securely.

Powder-free and Latex-free

Eliminates powder-induced irritation and makes it ideal for people with a Type I latex sensitivity.

Minimises the Risk of Allergic Contact Dermatitis

Free of Thiurams, Mercaptobenzotiazoles, Thiazoles and Thioureas.

Chemical Resistant

Tested in accordance with EN374-1:2016 for resistance to chemical substances. Tested against chemotherapy drugs according to ASTM D6978-05.



CAT III



2777

EN ISO 374-1: 2016/Type B



KPT

EN ISO 374-5: 2016



VIRUS

EN455



THIS IS A MEDICAL DEVICE



EU Compliance

- MDR (EU) 2017/745 Class I
- PPE Regulation (EU) 2016/425 CAT III
- EN455
- Food Safe Regulation (EC) No 1935:2004 and EC 10/2011

Recommended For

- Accident and Emergency
- Bio-science
- Dental
- Laboratory
- Pharmacy
- Wards

Nitrile Examination Gloves

TECHNICAL

- **Type:** Non-sterile, powder-free single use examination gloves
- **Material:** Nitrile Butadiene Rubber (NBR)
- **Colour:** Light Blue
- **Glove Shape:** Flat with a beaded cuff
- **Impermeability:** AQL 1.5
(AS PER EN 455-1)

- **Force at Break (N):** >6 N
(AS PER EN 455-2)
- **Shelf Life:** 5 years
(IN ORIGINAL PACKAGING WHEN STORED AS PER DIN 7716, ISO 2230)
- **Country of Origin:** China

Material Information

- **Premium Performance Nitrile Butadiene Rubber**
- **Latex free**
- **Powder-free**
- **Beaded cuff**
- **Light Blue**

Sizing

- **Overall Length Size:** 240mm
(AS PER EN455-2)
- **Wall Thickness:** Min. 0.07 mm
(MEASURED IN SINGLE LAYER OF THE PALM AREA)
- **Weight:** 3.5 gr (size M)

Ordering Information

- **Product Code:**

HAK5100	XS
HAK5101	S
HAK5102	M
HAK5103	L
HAK5104	XL
HAK5105	XXL
- **Pack Quantity:** 100 single gloves by weight, in dispensers
- **Case Quantity:** 10 dispensers per case

General Information

Conformance:

HAIKA™ NX510 medical examination gloves have been tested and found to comply with all parts of EN455.

Registering Authority: In Europe the gloves are CE marked indicating compliance with the PPE regulation 2016/425 (CAT III).

EU Type examination carried out by:
SATRA Technology Europe, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland. (Notified Body No. 2777)

They are also in compliance with Medical Device Regulation 2017/745 Medical Device Class 1 (non-sterile).

More information:

EUROPE

www.globusgroup.com

E: sales@globus.co.uk
T: +44 (0)161 877 4747
F: +44 (0)161 877 4746

MIDDLE EAST AND AFRICA

www.globusgroup.com/gcc

E: gcc@globusgroup.com
T: +971 4 882 9962
F: +971 4 882 9963

