

2- door tall safety cabinet Type 90 - 4 compartments

COMPLIANCE

EN 14470-1, EN 1363-1

ADVANTAGES

- 105-minute fire resistant and made in France
- EN 14470-1 and FM 6050 compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Thermo-expanding seals of the doors
- Locking system
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Fixation for earthing
- 90-minute fire resistant according to EN 14470-1, and tested until 105 minute
- Tested and approved according to FM 6050
- Tested and approved with ISO 834 fire heating curve by an European laboratory
- Compartments allowing the separation of different types of chemicals



TECHNICAL FEATURES

| | | | |
|-------------------|--|--------------------------------|--|
| Model | 795+M22 | Shipping dimensions HxWxD (mm) | 2150 x 1200 x 800 |
| Compliance | EN 14470-1, EN 1363-1 | Ventilation | Both inlet and outlet on the back wall of the cabinet |
| Cabinet materials | 12/10th epoxy painted steel structure | Openings diameter (mm) | 100 |
| Doors color | white RAL 9010 | Number of compartments | 4 |
| Number of doors | 2 | Containment sump | 4 |
| Doors type | swing doors | Grating | 4 |
| Hinges | right & left | | 497.00 |
| Door closure | self-closing | HS Code | 94032080 |

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Containment sump for 1-door cabinets Types 30, 60 & 90 (Réf : B148)
- Perforated cover plate for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : C148)
- Retention shelf for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : E48)
- Containment sump in PVC for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : SB148)
- Protection shelf in PVC for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : SE148)
- Sliding drawer for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : T148)

| Model | External dimensions HxWxD (mm) | Internal dimensions HxWxD (mm) | Storage capacity (L) | Containment sump | Shelves | Load capacity | Weight (kg) |
|---------|--------------------------------|--------------------------------|----------------------|------------------|---------|---------------|-------------|
| 795+M22 | 1950 x 1137 x 620 | (2 x 800) x (2 x 490) x 410 | 220 | | 4 | 30 | 497.00 |