

**RANGE 8B**

## 2- DOOR WORKING COVER STYLE MELAMINE SAFETY CABINET



### CONFORMITY

-

### ADVANTAGES

- Especially designed for acid and bases storage
- Built in melamine-based wood with a great resistance to corrosive vapours aggressions.
- Cabinets with sliding drawers for acid or/and base storage

### TECHNICAL CHARACTERISTICS

- Construction of the cabinet: Melamine wood structure
- Body color: white RAL 9010
- Door(s) color: white RAL 9010
- Number of door(s): 2
- Door styles: sliding doors
- Hinges:
- Fermeture de(s) porte(s) : manual closing
- Key locking system: manual closing
- Ventilation: Both inlet and outlet
- Outlet diameter (mm): 100

- Connecting sleeve included: included
- Number of compartment(s) : 2
- Shelf(ves) : -
- PVC containment shelf: -
- Containment sump: -
- PVC containment sump: -
- Perforated cover plate: -
- Shipping dimensions H x W x D (mm): 1100 x 1200 x 800
- Volume (m3) : 0,6282225
- Customs code : 94036000

### FILTERS

- Ventilation box (Réf : CDV-A)
- Filtration box with active coal filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box in PVC with active coal filter for organic and corrosive vapors (Réf : CDF-P)
- Recycling air box Labopur (Réf : H50C)

### OPTIONS

- Sliding drawer-containment sump (Réf : E8A2)

P/N	External dim. H x W x D (mm)	Internal dim. H x W x D (mm)	Storage capacity (L)	Retention capacity (L)	Number of shelves	Maximum load by shelf (kg)	Weight (kg)
APC112	900 x 1135 x 615	2x (760 x 465 x 535)	110	66	-	20	

**ECOSAFE**  
www.ecosafesa.com

Chemin des Champs Courbes - Zone Industrielle de  
Vallaire Est  
CH-1204 - Ecublens (Switzerland)  
Tel. +41 (0)21 694 10 04 - Fax : +41 (0)21 694  
10 00  
www.ecosafesa.com - sales@ecosafesa.com

Flammable  
cabinetsCorrosive  
cabinetsToxic  
cabinets  
and  
multiriskPesticide  
rangeFume  
hoods,  
filtering  
cabinets -  
VentilationContainment  
and  
CansAnti-fire  
equipments  
and Files  
cabinetsShowers  
and first-  
aid  
equipments

