

7.5L QUIVER

Product Code:

5190XXXX

Our Sharpsafe containers are specifically designed to enable the safe disposal of clinical sharps and minimise the risk of associated injuries.

Our 7.5L Quiver container is designed with operating theatres in mind to ensure the safe disposal of longer items used during procedures, such as endoscopy instruments. The container is also perfectly sized to store catheters.

Available with different lids for easily identifiable waste segregation and a number of compatible accessories for increased stability, security and safety.

Sharpsafe's new 5th Generation containers launched in 2022 and were the world's first full range of sustainable single use sharps containers. Engineered using ultra low carbon footprint recycled material showing our commitment to championing change.



SUPPORTING
A CIRCULAR
ECONOMY



REDUCING
YOUR CARBON
FOOTPRINT



REDUCING VIRGIN
PLASTIC ENTERING
THE SUPPLY CHAIN



MADE USING LOCALLY
SOURCED LOW
CARBON FOOTPRINT
RECYCLED MATERIALS



Product Certification	
ISO 23907-1:2019	✓
UN3291	✓
BSI Kitemark	✓

Site Certification	
ISO 9001 - Quality Management	✓
ISO 14001 - Environmental Management	✓
ISO 45001 - Occupational Health & Safety	✓

Product Details	
Size:	7.5L
Usable Capacity:	6.15L
Usable Volume:	84%
Weight:	880g/954g [^]
Fill Line Height:	500mm
Autoclavable:	Yes
Autoclavable Conditions:	134°C
Product Dimensions (mm):	610 x 150 x 100 (HxWxD)
Aperture Dimensions (mm):	140 x 90 (LxW)
GTIN:	5056271506111
Lot/Batch Number:	Not Applicable
Expiry Date:	Not Applicable

Accessory Compatibility			
Wall Bracket	n/a	Pedestal Trolley	n/a
Stand	✓	Sharpsafe Trolley	n/a
NPSD Tray	n/a	XL Trolley	n/a

Logistic Data	
Carton Quantity:	20
Carton Weight (when empty):	1.56kg
Carton Weight (when full):	19.5kg
Carton Dimensions (mm):	615 x 450 x 765 (HxWxD)
Carton Volume:	0.212m ³
Storage Temperature Conditions:	Ambient
Storage Lower Temp Limit:	None
Storage Upper Temp Limit:	None
Sterile/Non Sterile:	Non Sterile

Features and Benefits			
Self-Adhesive Label	✓	Latex Free	✓
Temporary Closure	✓	Low Carbon Footprint Material	✓
Permanent Closure	✓	Leak Resistance	✓
Four Click Assembly	n/a	Disconnection Features	n/a
Fill Line	✓	Safety Shield	✓
Protected Access*	✓	Semi-Translucent Lid	n/a
Single Use	✓	Carry Handle	✓
Different Lid Variations*	✓	Puncture Resistant	✓

Pallet Information	
Euro Pallet Carton Quantity:	6
Euro Pallet Dimensions (mm):	1200 x 800
Standard Pallet Carton Quantity:	9
Standard Pallet Dimensions (mm):	1200 x 1000

Manufacturing Information	
Made In:	England, UK
By:	Vernacare

[^]Weight without protected access insert/Weight including protected access insert
^{*}Different product codes may apply

FOR MORE INFORMATION:

✉ sharpsafesales@vernagroup.com

🌐 vernacare.com 📞 +44 (0)1495 235800