

Safety Unlimited  
420000  
Unit 2  
528128  
40 Comet Way  
Eastwoodbury Lane  
Southend  
Essex  
SS2 6XW

Tel: 01702

Fax: 01702

e-mail: [sales@safetyunlimited.co.uk](mailto:sales@safetyunlimited.co.uk)  
[www.safetyunlimited.co.uk](http://www.safetyunlimited.co.uk)

**Product Code:** SAF8090

**Description:** Radioactive storage cabinet for use with low level radioactive materials

**Safety  
Unlimited**

### Information:

Manufactured from 1.2mm Zintec coated sheet steel for excellent corrosion resistance.

Spot welded seams with gas welded corner joints.

Dimensions: 380mm x 380mm x 380mm (H x W x D)

One fixed, half depth shelf and plastic contents spray.

Labelled for easy identification.

Secure locking mechanism and louvered door.

Finished in grey, polyester powder coat.

Approximate weight: 8kg



### Use

This storage unit is designed to prevent environmental pollution and provide secure storage for low level radioactive substance.