


1. Part Number: 1210K-2000

2. Photo:



3. Product Information:

Description	Jumbo Slot Sampler 2000mm
Material of Construction:	316L stainless steel
Outer Diameter of Outer Tube:	40 mm
Inner Diameter of Inner Tube:	32 mm
Overall Length:	2060 mm
Slot Length:	120 mm
Slot Width:	20 mm
Number of Slots	4
Slot Positions:	



External Polish:	<1 Micron Ra
Internal Polish:	None
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Nominal Weight of Sampler:	6.5 Kg
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Oct-19