


**1. Part Number:** 8060H-500.5S

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare PowderThief 500mm – 0.5ml Tip - Sterile
Use:	Single use, disposable powder sampler
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Ultrasonic welding
Nominal Length:	511mm
Diameter of Outer Tube:	21mm
Nominal Tip Size:	0.5ml
Width of Handle:	103mm
Nominal Weight of Sampler:	93g
Assembly & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (minimum dose 25 kGy)
Number of Samplers per Box:	20
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Jun-20