

**1. Part Number:** 8085A-100

**2. Photograph:**



### 3. Product Information:

Description:	Disposable Sack Sampler – 400mm
Use:	Single use, disposable powder sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> </ul>
Method of Construction:	Single piece construction
Overall Length:	400mm
Diameter of Outer Tube:	25mm
Nominal Weight of Sampler:	53g
Quantity per Box:	100
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	N/A
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	2 years from the date of manufacture

### 4. Change Control:

Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

### 5. Country of Origin:

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21