


1. Part Number: 8073J-1000

2. Photograph:



3. Product Information:

| | |
|----------------------------|---|
| Description: | SteriWare Jumbo Slot Sampler - 1000mm |
| Use: | Single use, disposable powder sampler |
| Material of Construction: | High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free |
| Method of Construction: | Ultrasonic welding |
| Overall Length: | 986mm |
| Diameter of Outer Tube: | 40mm |
| Diameter of Inner Tube: | 35mm |
| Nominal Weight of Sampler: | 297g |
| Number of Sampling Slots: | 3 |
| Sampling Slot Width: | 8mm |
| Sampling Slot Length: | 60mm |
| Sampling Slot Positions: | |



| | |
|---|------------------------------|
| Assembly & Primary Packing Environment: | Class 100,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | N/A |



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

Number of Samplers per Box: 20
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21